


SysLoc[®] MINI

JMS SysLoc[®] Mini with Rotating Bevel


The latest design in safety AVF Needles. The SysLoc[®] MINI combines innovative technology and the audible "click" to ensure safety for patients and health-care workers.



SysLoc[®] MINI is over 20% shorter than the previous model providing easier use.



Shorter transparent hub improves control while providing full view of flash-back for proper positioning.



Contoured rotating hub provides for optimal flow and patient comfort.




The audible click ensures that the needle is securely locked into place – reducing the chance for needle-stick injury.



www.jmsna.net

JMS

System Safety Activation



For additional operating instructions please refer to the instructions for use.

Ordering Information

Single Set

Item Code	Description	Each / Case
864-1400-33	AVF Set 14G X 1" BE w/Clamp 12" Tubing	300
864-1500-33	AVF Set 15G X 1" BE w/Clamp 12" Tubing	300
864-1600-33	AVF Set 16G X 1" BE w/Clamp 12" Tubing	300
864-1700-33	AVF Set 17G X 1" BE w/Clamp 12" Tubing	300
864-1800-33	AVF Set 18G X 1" BE w/Clamp 12" Tubing	300
864-1504	AVF Set 15G X 1" 1/4 BE w/Clamp 12" Tubing	300
864-1605	AVF Set 16G X 1" 1/4 BE w/Clamp 12" Tubing	300
864-1704	AVF Set 17G X 1" 1/4 BE w/Clamp 12" Tubing	300

Twin Set

Item Code	Description	Twin Set / Case
864-2400-33	AVF Twin Set 14G X 1" BE w/Clamp 12" Tubing	250
864-2500-33	AVF Twin Set 15G X 1" BE w/Clamp 12" Tubing	250
864-2600-33	AVF Twin Set 16G X 1" BE w/Clamp 12" Tubing	250
864-2700-33	AVF Twin Set 17G X 1" BE w/Clamp 12" Tubing	250
864-2800-33	AVF Twin Set 18G X 1" BE w/Clamp 12" Tubing	250
864-2502	AVF Twin Set 15G X 1" 1/4 BE w/Clamp 12" Tubing	250
864-2602	AVF Twin Set 16G X 1" 1/4 BE w/Clamp 12" Tubing	250
864-2702	AVF Twin Set 17G X 1" 1/4 BE w/Clamp 12" Tubing	250

Imported by:

JMS North America Corporation

22320 Foothill Blvd., Suite 350
Hayward, CA 94541
Tel: 510-888-9090 • Fax: 510-888-9099

www.jmsna.net info@jmsna.net



[www.Facebook.com/JMSNA](https://www.facebook.com/JMSNA)



Distributor